



monolithic ZR with zero cutback

miYO<sup>®</sup>  
esthetic system

miYO<sup>®</sup>  
pink



Importers & Dealers in Quality Dental Goods

watch MiYO restoration  
being made here!  
[jensendental.com/see-miyo](http://jensendental.com/see-miyo)

patent pending

restorations by Don Cornell & James Choi  
photography by James Choi



Raise value without adding opacity

## MiYO Color

Translucent, semi-translucent and opaque self-glazing colors that don't run or puddle



sintered zirconia



MiYO Color wet



MiYO Color fired

Patent pending materials allow users to float color on top of wet color without mixing, so that every detail imaginable can be achieved in one application



before firing



MiYO Structure Window was used to detail shape and texture (before firing)

## MiYO Structure

Structure building materials in varying colors and translucency for the purpose of building shapes or detailing surface texture, for teeth and tissue



after firing



All details faithfully preserved and fully glazed



# No Cutback

Achieve depth, translucency and effects in 0.1 – 0.2mm with MiYO Color & Structure



**miYO**  
reimagine monolithic.  
visit [jensendental.com/see-miyo](http://jensendental.com/see-miyo)

**JENSEN**  
**Chemichl**  
Custom Designed Ceramics  
Manufactured by Chemichl AG, an affiliate of Jensen Dental in Vöcklabruck, Liechtenstein



Importers & Dealers in Quality Dental Goods

1102 Embassy Chambers, 3rd road, Khar West, Mumbai 400052 
 +91 22 6127 7777  
[sales@nkpatelandsons.com](mailto:sales@nkpatelandsons.com)
 +91 99300 06127  
[@nkpatelandsons.com](https://www.facebook.com/nkpatelandsons.com)
[/nkpatelandsons.com](https://www.instagram.com/nkpatelandsons.com)